

## MTG Adolescent & Kiddie-Kath



ISC closed system designed especially for children.

- Pre-lubricated 10" catheter
- Self contained 700ml collection bag
- Available as a discreet Mini-Pak
- Privacy disposal bag included

Description	NHS Order Code	Size	PIP Code
MTG EZ-Advancer Adolescent	30110	10	395-0516
MTG Kiddie Kath	30108	8	395-0508

## MTG EZ-Gripper

Simply the easiest closed system ISC for patients with compromised hand function or limited dexterity. Revolutionary patented gripper device that aids insertion of the catheter. A real alternative to an indwelling catheter.



- Pre-lubricated 16" catheter
- Self contained 1500ml collection bag with attached outer gripping device
- Privacy disposal bag
- The only catheter that can be used with the MTG Eagle Board

Description	NHS Order Code	Size	PIP Code
MTG EZ-Gripper	50112	12	
MTG EZ-Gripper	50114	14	
MTG EZ-Gripper	50116	16	



# MTG RANGE OF ISC CATHETERS

Effortless worry-free catheterisation to fit your lifestyle



<https://www.facebook.com/PioneerMedicalEurope>



<https://twitter.com/PioneerEU>

Please phone for samples 0800 0831 337 or email [samples@mtgcatheters.com](mailto:samples@mtgcatheters.com)

[www.mtgcatheters.com](http://www.mtgcatheters.com)

**PIONEER**

Brought to you by  
Pioneer Medical Europe distributor  
for the UK and Ireland

Unit 340 Beech Drive, Hartlebury Trading Estate,  
Hartlebury, Worcestershire, DY10 4JB Tel: 01299 489847

UK Freephone: 0800 0831 337 International: 0044 (0) 845 303 5868

email [sales@pioneermedicaleurope.co.uk](mailto:sales@pioneermedicaleurope.co.uk)

Pioneer Medical Europe Ltd [www.pioneermedicaleurope.co.uk](http://www.pioneermedicaleurope.co.uk)

**PIONEER**

Brought to you by  
Pioneer Medical Europe distributor for  
the UK and Ireland



# MTG Range of Catheters

Protect the sterility of the catheter during the entire procedure and collect the urine upon voiding

All MTG catheters offer:

- MTG soft silicone Introducer Tip is designed to shield the catheter from contact with bacteria located within the first 1.5cm (½ inch) of the distal urethra. So if the catheter is passed beyond this region, without coming into contact with these potential pathogens, UTI recurrence risk is reduced
- The patent pending MTG EZ-Advancer valve is a non-return mechanism that secures the catheter tube after each forward advancement, preventing retraction during insertion and tube slippage back into the bag. The unique locking mechanism enables very easy insertion
- MTG range is easy to use. Unlike other ISC catheters, these catheters are all designed especially for users with dexterity impairment
- MTG catheter range is a closed system, allowing touch-free catheterisation. The bag serves to protect the sterility of the catheter and collect fluid
- Easy open packaging
- Fire Polished eyelets for comfort
- 100% Latex Free & DEHP Free
- Available male, female and childrens

In summary, MTG closed system is intended to protect the sterility of the catheter during the entire procedure and collect the urine upon voiding.

## MTG EZ-Advancer



Available as a standard pack or Mini-Pak for those on the move.

- Pre-lubricated 16" catheter
- Self contained 1500ml collection bag
- Privacy disposal bag included in Mini-Pak versions (40112, 40114, and 40116)

Description	NHS Order Code	Size	PIP Code
MTG EZ-Advancer (Male Soft)	30212	12	395-0409
MTG EZ-Advancer (Male Soft)	30214	14	395-0417
MTG EZ-Advancer (Male Firm)	30112	12	395-0425
MTG EZ-Advancer (Male Firm)	30114	14	395-0433
MTG EZ-Advancer (Male Firm)	30116	16	395-0441
MTG EZ-Advancer Mini Pak (Male)	40112	12	395-0458
MTG EZ-Advancer Mini Pak (Male)	40114	14	395-0466
MTG EZ-Advancer Mini Pak (Male)	40116	16	395-0474

## MTG Cath-Lean



ISC closed system designed especially for women.

- Discreet and inconspicuous
- Pre-lubricated 6" length female catheter
- Self contained 700ml collection bag
- Available as a Mini-Pak, fits easily in handbag
- Privacy disposal bag included

Description	NHS Order Code	Size	PIP Code
MTG Cath-Lean (Female)	40412	12	395-0482
MTG Cath-Lean (Female)	40414	14	395-0490

Please phone for samples 0800 0831 337 or  
email [samples@mtgcatheters.com](mailto:samples@mtgcatheters.com)

[www.mtgcatheters.com](http://www.mtgcatheters.com)

Please phone for samples 0800 0831 337 or  
email [samples@mtgcatheters.com](mailto:samples@mtgcatheters.com)

[www.mtgcatheters.com](http://www.mtgcatheters.com)



**PIONEER**  
Brought to you by  
Pioneer Medical Europe distributor for  
the UK and Ireland